ACCERO® Rex Stent

ENGINEERING STROKE SOLUTIONS

Acandis GmbH

Theodor-Fahrner-Str. 6
75177 Pforzheim/Germany
www.acandis.com

DE Patienteninformation
EN Patient information

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1 Product Information

The ACCERO® Rex Stent is a self-expanding implant made of nitinol. In the centre, the implant has a radiopaque marker made of platinum-iridium.

2 Intended Use

2.1 Intended Purpose

The ACCERO® Rex Stent is used to anchor the embolisation materials in the aneurysm when using the stent-assisted coiling method.

2.2 Indication

The ACCERO® Rex Stent is intended for use with embolisation materials in the treatment of intracranial aneurysms.

3 Information on Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated that the ACCERO® Rex Stent is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial magnetic field gradient of 7.0 T/m at 3.0 T and 11.0 T/m at 1.5 T
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) for 15 minutes of scanning 1.7 W/kg at 1.5 T and 2.7 W/kg at 3 T

Under the scan conditions defined above, the ACCERO® Rex Stent is expected to produce a maximum temperature rise of less than 3.9 °C at 1.5 T/2.8 °C at 3.0 T after 15 minutes of continuous scanning.

In non-clinical testing, the image artefact caused by the implant extends up to 8 mm from the ACCERO® Rex Stent when imaged with a spin echo pulse sequence and a 3.0 T MR system. Optimisation of the imaging parameters is recommended. In non-clinical testing, the magnetically induced displacement force and magnetically induced torque were tested and no clinically significant displacement or torque was measured.

The MR compatibility of the implant has only been established with a single stent. Acandis® offers no guarantee as to the MR compatibility of overlapping stents.

Please bring this information together with your implant card to all future MRI appointments.

4 Further Information

The ACCERO® Rex Stent is implanted for life. The ACCERO® Rex Stent is a passive implant that does not require maintenance or monitoring. Nevertheless, it is necessary for you to follow the instructions of your treating physician. These include taking the prescribed medications and attending arranged follow-up examinations. In combination with a healthy lifestyle, this can reduce the risk of recurring symptoms.

If you have any questions, medical concerns or health-related symptoms and with regard to recommended action and required follow-up examinations, please consult your treating physician or your specialist care team.

In an emergency, please notify the emergency medical service.

Your treating physician will inform you when you should seek medical assistance. Furthermore, he will let you know of any potential side effects that might occur in connection with the implant, such as: unconsciousness, headaches, drowsiness, dizziness, speech disorders, uneven pupil dilation, nausea, fluctuations in blood pressure, mental confusion. Please discuss with your treating physician which of these side effects might occur in your case.

This implant was manufactured in consideration of all requirements relating to medical devices. Biocompatibility tests have shown that any manufacturing residuals do not compromise the safety of the product.

5 Reporting of Adverse Events

Please consult your treating physician to report side effects and adverse events that might be related to the implant.

For more information on this implant and updates to the information on this printout, please visit our website at: www.acandis.com/productinfo

You can report potential serious events in connection with this implant to Acandis®. Additionally, patients in Australia can contact the Australian Therapeutic Goods Administration (TGA) at: www.tga.gov.au